

HomeFill® Cylinders

M2 Cylinder

ML4 Cylinder

ML6 Cylinder

M9 Cylinder

Integrated Conserver Cylinders



Model Number	HF2PCM2KIT	HF2PCL4	HF2PCL6	HF2PC9
Height	8.25"	9"	11"	14"
Weight	2 lb	3.1 lb	3.5 lb	4.8 lb
Diameter	3"	4.38"	4.38"	4.38"
Capacity	57 Liters	113 Liters	170 Liters	255 Liters

Conserved Durations (Hrs)*				
At setting 1	3.4 Hrs	6.7 Hrs	9.8 Hrs	15.2 Hrs
At setting 2	2.2 Hrs	4.3 Hrs	6.3 Hrs	10.5 Hrs
At setting 3	1.6 Hrs	3.2 Hrs	4.7 Hrs	7.3 Hrs
At setting 4	1.3 Hrs	2.6 Hrs	3.8 Hrs	5.9 Hrs
At setting 5	1.1 Hrs	2.1 Hrs	3.1 Hrs	4.8 Hrs
Fill Time	24 Min	55 Min	75 Min	125 Min
Flow Settings	1-5 Pulse or 2 lpm Cont.	1-5 Pulse or 2 lpm Cont.	1-5 Pulse or 2 lpm Cont.	1-5 Pulse or 2 lpm Cont.

*20 bpm

Post Valve ML6

Post Valve M9

Post Valve D

Post Valve Cylinders



Model Number	HF2POST6	HF2POST9	HF2POSTD
Height	12"	15"	21"
Weight	4.2 lb	5.4 lb	6.5 lb
Diameter	4.38"	4.38"	4.38"
Capacity	170 Liters	255 Liters	422 Liters
Duration @ 2 lpm flow	Varies	Varies	Varies
Fill Time	75 Min	125 Min	3.5 Hrs
Flow Settings	N/A	N/A	N/A

M6 Cylinder

M9 Cylinder

D-Tank

Continuous Flow Cylinders



Model Number	IOH2M6	IOH2M9	IOH2MD
Height	15.25"	14.25"	19.75"
Weight	3.5 lb	4.8 lb	7 lb
Diameter	3"	4.38"	4.38"
Capacity	148 Liters	248 Liters	422 Liters
Duration @ 2 lpm flow	1 Hr	1.9 Hrs	3.5 Hrs
Fill Time	1 Hr	1.9 Hrs	3.5 Hrs
Flow Settings	0 - 6 lpm	0 - 6 lpm	0 - 6 lpm

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44036-4190
(800) 333-6900

Canada
570 Matheson Blvd E Unit 8
Mississauga Ontario
L4Z4G4 Canada
(800) 668-5324

©2010 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Form No. 07-066 Rev.08/10 100433

